


eMorpho-HPS 2-inch square embedded MCA



2"x2" square PC-board
Origin at lower left hand corner of PCB

Connector	Center (x,y)	Extent (x,y)
USB	(0.105", 0.260")	(0.356", 0.290")
GPIO	(0.590", 0.670")	(0.866", 0.197")
Detector	(0.275", 1.220")	(0.197", 0.551")
JTAG	(0.360", 1.830")	(0.655", 0.246")

Mounting holes: 0.100" diameter

#1	(0.100", 1.550")
#2	(1.825", 1.850")
#3	(1.850", 0.150")

eMorpho-HPS outline drawing